



Specifications

Product : SOYABEAN OIL
CAS No. : 8001-22-7
Product Code : 02NA0141

TESTS	SPECIFICATIONS
Color and appearance	A pale yellow oil
Odor	Odorless or almost odorless.
Assay	It is extracted from the seeds of the soybean (Glycine max) using n-hexane as solvent contains 22% of the dry seed.
Solubility	Practically insoluble in ethanol (96%); miscible with chloroform, with ether and with petroleum spirit
Acid value	NMT 0.6
Refractive index@20°C	1.465 to 1.475
Wt. per ml@25°C	0.916 to 0.922
Unsaponifiable Matter	NMT 1.5%
Stearin	Remains bright when maintained at 0° for 16 hours.
Foreign fixed oil	Palmitic acid: 7.0% to 14.0% Stearic acid: 1.4% to 5.5% Oleic acid: 19.0% to 30.0% Linoleic acid: 44.0% to 62.0% Linolenic acid: 4.0% to 11.0%

NMT= Not more than, NLT=Not less than)

INFORMATION	PREPARED BY	APPROVED BY
Name	Mr. Suchit Mukherjee	Mr. Thomas Varghese
Designation	Q.C Executive	Director Technical
Date	Jan 2019	-
Revision	0	-
Review Date	-	-
Signature		

Disclaimer: The products sold hereunder and any other products sold by Indenta Chemicals (India) Pvt. Ltd. have been subjected to limited testing and it does not release the user from the responsibility of undertaking their own test.

Indenta Chemicals (India) Pvt. Ltd.

Office: 117 The Summit Business Bay, Near WEH Metro Station, Opp. Cinemax Theatre, Off. Andheri Kurla Road, Andheri (E), Mumbai 400 093.
 Phone : +91-22-2684 9600 | Fax : +91-22-2684 9060 | Email: indenta@indenta.com | Website : www.indenta.com

Unit 1: Plot No. 1405, GIDC Sarigam, Dist. Valsad, Gujarat – 396155

Unit 2: Building No. 73, Gala No. 7, Indian Corporation Compound, Village Gundavli, Mankoli Naka, Bhiwandi, Thane - 421302