



Specifications

Product : LEMON OIL
CAS No. : 8008-56-8
Product Code : 03GS0123

TESTS	SPECIFICATIONS
Color and appearance	A clear, mobile, pale yellow to dark green liquid, free from sediment & suspended matter.
Assay	Citral (California- type lemon oil): 2.2% - 3.8% Citral (Italian- type lemon oil: 3.0% -5.5%
Specific Gravity@25°	0.849 – 0.855
Refractive Index@20°	1.473 – 1.476
Optical Rotation	+57.0°to + 65.6°
Ultraviolet absorbance	The absorbance, calculated on the basis of a 250-mg specimen, is NLT 0.20 for California type lemon oil or NLT 0.49 for Italian type lemon oil.
Foreign oils	The angular roatation of distilled oil is NMT 6° less than that of the original oil.

(NMT= Not more than, NLT=Not less than)

INFORMATION	PREPARED BY	APPROVED BY
Name	Mr. Suchit Mukherjee	Mr. Thomas Varghese
Designation	Q.C Executive	Director Technical
Date	Jan 2019	-
Revision	0	-
Review Date	-	-
Signature		

Disclaimer: The products sold hereunder and any other products sold by Indenta Chemicals (India) Pvt. Ltd. have been subjected to limited testing and it does not release the user from the responsibility of undertaking their own test.

Indenta Chemicals (India) Pvt. Ltd.

Office: 117 The Summit Business Bay, Near WEH Metro Station, Opp. Cinemax Theatre, Off. Andheri Kurla Road, Andheri (E), Mumbai 400 093.
 Phone : +91-22-2684 9600 | Fax : +91-22-2684 9060 | Email: indenta@indenta.com | Website : www.indenta.com

Unit 1: Plot No. 1405, GIDC Sarigam, Dist. Valsad, Gujarat – 396155

Unit 2: Building No. 73, Gala No. 7, Indian Corporation Compound, Village Gundavli, Mankoli Naka, Bhiwandi, Thane - 421302

