

Product: LABDANUM RESINOID PURE-200142 **Code:** 200142

Chemical Name:

General product information

Natural Product.

CAS: 84775-64-4

EINECS: 283-893-1

CE Index-Nr.:

Technical specification

TEST	LIMIT VALUES	UNITS	METHOD
Appearance:	Paste-semi-solid	---	QC 03
Colour:	Dark brown	---	QC 03
Solubility at 20 °C in alcohol 94%:	Turbid	---	QC 04
Flash Point:	+100	°C	(TCC) QC 08

Additional information

Odour: Characteristic.

Storage: In full tight sealed containers at temperature below +20 °C and protected from light.

Shelf Life: 24 months.