



Specifications

Product	:	Methyl salicylate
CAS No.	:	119-36-8
Molecular formula	:	C ₈ H ₈ O ₃
Molecular weight	:	152.149
Product Code	:	MC0043

TESTS	SPECIFICATIONS
Color and appearance	A clear colorless liquid.
Assay	99.0-100.5%
Appearance of solution	The solution is clear
Acidity	NMT 0.4ml of 0.1M NaOH is required to restore the blue color.
Refractive index	1.535-1.538
Specific gravity	1.180-1.186
Heavy metal	NMT 20ppm

(NMT= Not more than, NLT=Not less than)

INFORMATION	PREPARED BY	APPROVED BY
Name	Mr. Suchit Mukherjee	Mr. Thomas Varghese
Designation	Q.C Executive	Director Technical
Date	Jan 2019	-
Revision	0	-
Review Date	-	-
Signature		

Disclaimer: The products sold hereunder and any other products sold by Indenta Chemicals (India) Pvt. Ltd. have been subjected to limited testing and it does not release the user from the responsibility of undertaking their own test.

Indenta Chemicals (India) Pvt. Ltd.

Office: 117 The Summit Business Bay, Near WEH Metro Station, Opp. Cinemax Theatre, Off. Andheri Kurla Road, Andheri (E), Mumbai 400 093.
 Phone : +91-22-2684 9600 | Fax : +91-22-2684 9060 | Email: indenta@indenta.com | Website : www.indenta.com

Unit 1: Plot No. 1405, GIDC Sarigam, Dist. Valsad, Gujarat – 396155

Unit 2: Building No. 73, Gala No. 7, Indian Corporation Compound, Village Gundavli, Mankoli Naka, Bhiwandi, Thane - 421302

