



Specifications

Product	:	Polyethylene glycol 6000
CAS No.	:	25322-68-3
Molecular formula	:	$C_{2n}H_{4n+2}O_{n+1}$
Molecular weight	:	6000
Product Code	:	MC0045

TESTS	SPECIFICATIONS
Color and appearance	A white to creamy white, wax-like solid or flakes.
Appearance of solution	15.0% w/v solution is not more intensely colored than reference solution BSY6.
pH	4.5-7.5
Freezing point	56-60°C
Viscosity	250-390mm ² S ⁻¹
Arsenic	NMT 3 ppm
Heavy metals	NMT 5 ppm
Sulphated ash	NMT 0.1%

(NMT= Not more than, NLT=Not less than)

INFORMATION	PREPARED BY	APPROVED BY
Name	Mr. Suchit Mukherjee	Mr. Thomas Varghese
Designation	Q.C Executive	Director Technical
Date	Jan 2019	-
Revision	0	-
Review Date	-	-
Signature		

Disclaimer: The products sold hereunder and any other products sold by Indenta Chemicals (India) Pvt. Ltd. have been subjected to limited testing and it does not release the user from the responsibility of undertaking their own test.

Indenta Chemicals (India) Pvt. Ltd.

Office: 117 The Summit Business Bay, Near WEH Metro Station, Opp. Cinemax Theatre, Off. Andheri Kurla Road, Andheri (E), Mumbai 400 093.
 Phone : +91-22-2684 9600 | Fax : +91-22-2684 9060 | Email: indenta@indenta.com | Website : www.indenta.com

Unit 1: Plot No. 1405, GIDC Sarigam, Dist. Valsad, Gujarat – 396155

Unit 2: Building No. 73, Gala No. 7, Indian Corporation Compound, Village Gundavli, Mankoli Naka, Bhiwandi, Thane - 421302

